

100% Made in Germany

Products of the highest Standard



www.mk-medical.com



mk medical
orthopedic implants

Content / Delivery program

Kirschner Wires & Steinmann Pins.....	3	Dispenser for K-Wires	13
MK Numbering system for Kirschner Wires & Steinmann Pins..	5	Flexible Dispenser for K-Wires	15
The Tips in Detail.....	6	K-Wire Dispenser Set	16
Kirschner Wires.....	7	Storage Rack for Kirschner Wires	17
Titanium Kirschner Wires.....	7	Sterilization Containers	18
Steinmann Pins.....	7	Sterilization Containers for Kirschner Wires.....	19
Titanium Steinmann Pins	7	Caliper for Kirschner Wires	19
Bone Soft Wire	8		
Bone Soft Wire from Steel or Titanium	9		
Intramedullary Pins	10		
Intramedullary Pins made of Titanium	11		
Accessories	12		

Kirschner Wires & Steinmann Pins

PUTTING SAFETY FIRST





The latest generation of Kirschner Wires

Triple-S – precision machined with innovative, cutting edge technology.
Exclusively and directly from your leading OEM manufacturer.



Love it in color? No Problem!

Exclusive and only with us: Steel-Wires in your desired color.



MK Numbering system for Kirschner Wires & Steinmann Pins

Our article numbers build logically on to the following system:

T KB 16.150 TR.TG15

Item: Kirschner Wires (T=Titanium)

Diameter of the wire: 1,6 mm

Design: Tip 1/Tip 2. Special-Design

Length of the wire:150 mm



ITEM

KB Kirschner Wires

SN Steinmann Pin

Designs

T Trocar

R Round

F Flat

L Lancet

LB Lancet with hole

3 Triangle

4 Square

Special-Designs

TGx Part Thread
[x= Thread length in mm]

VG Full Thread

LR Trocar knurled lengthwise

OE Bolt Eye

B Drill Point

Custom designs available!

The listing is only a small extract from our manufacturing program. In case you cannot find the desired wire here please ask us. We shall be happy to give you detailed information.

The Tips in Detail

Standard-Tips



T Trocar



R Round



F Flat



L Lancet



LB Lancet with hole



3 Triangle



4 Square



B Drill Point


Special Designs



TG Partial Thread 15 mm



LR Trocar knurled lengthwise



VG Full Thread



OE Bolt Eye

Kirschner Wires

Ø 0,7 - 3,0 mm

Implant Steel 1.4441

Lengths of 50 - 600 mm

Packaging Unit = 10 pieces

Titanium Kirschner Wires

Ø 1,0 - 3,0 mm

Lengths from 50 - 600 mm

Generally all our Titan wires are available in all Standard Designs and with thread.

The Titan wires can be supplied in the following colours. Please specify the desired colour by indicating the number when placing your order.

Steinmann Pins

Ø 3,2 - 6,0 mm

Implant Steel 1.4441

Lengths of 50 - 600 mm

Packaging Unit = 1 piece

Titanium Steinmann Pins

Ø 3,5 - 5,0 mm

Lengths from 50 - 600 mm



Bone Soft Wire



Bone Soft Wire from Steel or Titanium

Ø 0,2 - 2,0 mm

Implant Steel 1.4441



Diameters	PU	Implant Steel	Titanium
Ø 0,2 mm	Coil, 10 m	MK-3020	
Ø 0,3 mm	Coil, 10 m	MK-3030	
Ø 0,4 mm	Coil, 10 m	MK-3040	
Ø 0,5 mm	Coil, 10 m	MK-3050	
Ø 0,6 mm	Coil, 10 m	MK-3060	
Ø 0,7 mm	Coil, 10 m	MK-3070	
Ø 0,8 mm	Coil, 10 m	MK-3080	TMK-3080
Ø 0,9 mm	Coil, 10 m	MK-3090	
Ø 1,0 mm	Coil, 10 m	MK-3100	TMK-3100
Ø 1,2 mm	Coil, 10 m	MK-3120	TMK-3120
Ø 1,5 mm	Coil, 10 m	MK-3150	TMK-3150
Ø 1,8 mm	Coil, 10 m	MK-3180	
Ø 2,0 mm	Coil, 10 m	MK-3200	

Intramedullary Pins



Intramedullary Pins made of Titanium

Length: 450 mm

Diameter	Item-ID				
Ø 1,5 mm	TIS15.450	3	6	9	12
Ø 2,0 mm	TIS20.450				
Ø 2,5 mm	TIS25.450	15	18	21	24
Ø 3,0 mm	TIS30.450				
Ø 3,5 mm	TIS35.450	27	30	33	36
Ø 4,0 mm	TIS40.450				
Ø 4,5 mm	TIS45.450	39	42	45	48
Ø 5,0 mm	TIS50.450	51	54	57	60
		63	66	69	72
		75	78	81	84
		87	90	93	95

The pins can be supplied in the adjacent colours. Please specify the desired colour by indicating the number when placing your order.

Other length 100 - 450 mm
and diameter 1,5 - 5,0 mm



Also available in stainless steel on request

Accessories



Dispenser for K-Wires

Ø 0,8- 1,4 mm, stainless steel


for optimum storage and cleaning of K-Wires up to a length of 150 mm.

Inspection hole for quick control of wire diameter, PTFE impact protector on both sides.

Classification: Ø 0,8/1,0/1,2/1,3/1,4 mm

Measurement: 162 x 74 x 40 mm

BS1164



Dispenser for K-Wires

Ø 1,4 - 3,0 mm, stainless steel


for optimum storage and cleaning of K-Wires up to a length of 150 mm.

Inspection hole for quick control of wire diameter, PTFE impact protector on both sides.

Classification: Ø 1,4/1,5/1,6/1,8/2,0/2,5/3,0 mm

Measurement: 162 x 96 x 40mm

BS1162



Dispenser for K-Wires

Ø 0,8 - 1,4 mm, stainless steel


for optimum storage and cleaning of K-Wires up to a length of 310 mm.

Inspection hole for quick control of wire diameter, PTFE impact protector on both sides.

Classification: Ø 0,8/1,0/1,2/1,3/1,4 mm

Measurement: 326 x 74 x 40 mm

BS1168



Dispenser for K-Wires

Ø 1,4 - 3,0 mm, stainless steel

for optimum storage and cleaning of K-Wires up to a length of 310 mm.

Inspection hole for quick control of wire diameter, PTFE impact protector on both sides.

Classification: Ø 1,4/1,5/1,6/1,8/2,0/2,5/3,0 mm

Measurement: 326 x 96 x 40 mm

BS1166



Flexible Dispenser for K-Wires

Ø 1,4 - 3,0 mm, stainless steel


for optimum storage and cleaning of K-Wires up to a length of 310 mm.

Inspection hole for quick control of wire diameter, PTFE impact protector on both sides.

Classification: Ø 1,4/1,5/1,6/1,8/2,0/2,5/3,0 mm

Measurement: flex. x 96 x 40 mm

BS1170



K-Wire Dispenser Set

Ø 0,8 - 3,0 mm; 150 mm

for optimum storage and sterile cleaning of K-Wires up to a length of 150 mm.

Integrated caliper for quick control of wire diameter,
PTFE impact protector on one side.

Classification: Ø 0,8/1,0/1,2/1,4/1,6/1,8/2,0/2,2/2,5/3,0 mm

Measurement: 236 x 74 x 38 mm

BS7600

Scope of Delivery:

1 x Sieve Tray


BS8016

1 x Lid

BS8015

1 x Wire Dispenser Insert
223 x 61 x 25 mm

BS1160



Storage Rack for Kirschner Wires


Ø1,0 - 3,0 mm

for optimum storage of K-Wires with a length of 150 mm.

Caliper holes for quick control of wire diameter,
silicon impact protector on both sides.

Ø1,0/1,2/1,4/1,6/1,8/2,0/2,5/3,0 mm

BS1500




for optimum storage of K-Wires with a length of 310 mm.

Caliper holes for quick control of wire diameter,
silicon impact protector on both sides.

Ø1,0/1,2/1,4/1,6/1,8/2,0/2,5/3,0 mm

BS3100



Sterilization Containers

NEW: Optimised liquid outlet



Sterilization Containers for Kirschner Wires

for wires from Item-ID.

100 mm **MK-1910**



160 mm **MK-1916**



230 mm **MK-1923**



310 mm **MK-1931**



440 mm **MK-1944**



other lengths available on request



with Tip Protection in Lid

Caliper for Kirschner Wires

for wires from Item-ID.

0,7 - 3,0 mm **MK-1950**





mk medical
orthopedic implants

MK medical orthopedic implants

Jahnstrasse 12

78576 Emmingen-Liptingen

Germany

Fon +49 7465 32601-0

Fax +49 7465 32601-99

info@mk-medical.de



www.mk-medical.com